

# Adaptive optics for microscopy

svmPSF quickstart manual  
vr2.0

doi: [10.5281/zenodo.3901390](https://doi.org/10.5281/zenodo.3901390)

[aomicroscopy.org](http://aomicroscopy.org)

Raphaël Turcotte, Carla C. Schmidt, Nigel J. Emptage, and Martin J. Booth  
Dynamic Optics and Photonics Group  
Department of Engineering, University of Oxford



## 1. DESCRIPTION

This user manual explains how to install and use the svmPSF plugins for Fiji/ImageJ. The svmPSF plugin models a spatially variant PSF into a series of spatially invariant PSFs. Details information about the svmPSF plugin and the algorithm it employs can be found in the [associated publication](#). Combined with a modified Richardson-Lucy deconvolution algorithm, the outputs from svmPSF enable deconvolution of images acquired with a multimode-optical fibre imaging system ([link to deconvolution code](#)).


If you use svmPSF for your research or education, please cite our [associated publication](#):

Turcotte R, Sutut E, Schmidt CC, Emptage NJ, Booth MJ.  
"Deconvolution for multimode fiber imaging: modeling of spatially variant PSF",  
Biomedical Optics Express 2020;11(8):759-4771. DOI: 10.1364/BOE.399983

## 2. INSTALLATION

To install the svmPSF plugin for Fiji/ImageJ follow these steps:

- 1) Click on the following link to access the svmPSF plugin: <https://github.com/dop-oxford/svmPSF>.
- 2) Download the svmPSF plugin file onto your computer (svmPSF\_1.0.0.jar).
- 3) Close any instance of Fiji/ImageJ.
- 4) Open the Fiji/ImageJ's plugin folder on your computer. Note: On MacOSX, use "Show Package Contents" to see the plugin folder.
- 5) Move the svmPSF plugin file into the Fiji/ImageJ's plugin folder.
- 6) Open Fiji/ImageJ. The svmPSF plugin should be listed under plugins (Fig. 1).



**Figure 1.** svmPSF plugin installation.


### 3. OPERATION

The operation of svmPSF can be tested with the provided datasets by following these steps:

- 1) Download and uncompress a test [dataset](#).
- 2) Open Fiji/ImageJ.
- 3) Click on *Plugins*, then *svmPSF* (Fig. 1).
- 4) Select the number of modes ( $N$ ).  $\sqrt{N}$  corresponds to the height and width of the PSF area. Use  $\sqrt{N} = 15$  with the test dataset.
- 5) Select the input folder within the dataset folder, which should include the point response measurements, and the output folder.
- 6) Select the position file. Note: It should be located in the input folder.
- 7) Click *Run*.

*Note: Running the svmPSF plugin will take approx. 30 seconds. The Run button stays blue as long as the plugin is running.*

- 8) Close the svmPSF plugin. The processed data is in the predefined output folder.



**Figure 2.** svmPSF plugin operation.